

Pull-down of RNA binding proteins with immobilized biotinylated transcripts

Buffers

- Cell lysis buffer:
 - 85 mM M KCl
 - 0.5% NP-40 (Pierce)
 - 5 mM HEPES, pH 7.4
 -
- Nuclei lysis buffer:
 - 300 mM NaCl
 - 20 mM Tris-HCl, pH 8.0
 - 0.2 % Tween-20
 - 1 mM EDTA
 - 1 mM EGTA
- 2x Binding/Washing buffer:
 - 10 mM Tris-HCl, pH 7.5
 - 1 mM EDTA
 - 2 M NaCl
 - 0.1 % Tween-20
- Streptavidin beads (Dynabeads)

1. Preparation of biotinylated transcripts

Materials

- Ambion MEGAscript T7 RNA-Polymerase Kit
- Bio-UTP (10 mM)
- Linearized template
- PCR machine
- RNeasy Kit (or Tri/TRIzol for purification of transcript).

Specific comment

- Original protocol from the Ambion Kit was changed to reduce the amount of biotinylated UTP.

Protocol

- Reaction mix (pipette in the order below):
 - 2 mM ATP, CTP and TTP

Protein pull-down with immobilized RNA

- 1 mM UTP
 - 1 mM bio-UTP
 - 2 µl 10x buffer
 - 1 µg linearized vector (use less in case of PCR products, see protocol from the kit).
 - 2 µl enzyme mix
 - H₂O ad 20 µl
- Incubate overnight at 37°C ??? (have to look up)
 - Add 1 µl DNase I from Ambion MEGAscript Kit, incubate at 37°C for 15 min.
 - Purify RNA with RNeasy Kit according to protocol supplied with the kit
 - Elute RNA with 50 µl H₂O

2. Immobilization of the biotinylated RNA to Streptavidin Dynabeads

Materials

- Streptavidin beads (Dynabeads)
- 2x Binding/Washing buffer:
 - 10 mM Tris-HCl, pH 7.5
 - 1 mM EDTA
 - 2 M NaCl
 - 0.1 % Tween-20

Protocol

Preparation of beads

- Wash Streptavidin beads (Dynabeads) 3x with Binding/Washing buffer (use 50 µl of bead suspension in a final assay with 1 ml of cell lysate)
- Wash beads 2x for 2 min with 50 µl of 0.1 M NaOH/0.05 M NaCl
- Wash once with the same volume of 0.1 M NaCl
- Resuspend beads in 50 µl of 0.1 M NaCl
- Transfer equal amounts of beads to DNA Low Binding tubes (0.5 ml sized)

Binding of RNA to beads

- Dilute biotinylated transcripts with H₂O (should have the same volume as the bead suspension)
- Denature RNA at 65°C for 10 min
- Cool down to RT over 20 min
- Add transcripts to beads; add water or unrelated RNA to control beads
- Incubate on a turning wheel at RT for 15 min
- Remove supernatant
- Measure OD₂₆₀ of the supernatant to verify efficient binding of the RNA to the beads
- Wash beads 3x with 1x Binding/Washing buffer

3. Preparation of cell lysates and pull down with immobilized RNA

Materials

- Mammalian cells growing adherently on 10 cm or 15 cm dishes
- Cell lysis buffer:
 - 85 mM M KCl
 - 0.5% NP-40 (Pierce)
 - 5 mM HEPES, pH 7.4
 -
- Nuclei lysis buffer:
 - 300 mM NaCl
 - 20 mM Tris-HCl, pH 8.0
 - 0.2 % Tween-20
 - 1 mM EDTA
 - 1 mM EGTA

Protocol

- *Extract preparation*
 - Remove medium, harvest cells with 1x PBS, transfer cells into a tube
 - Spin for 5 min at 800 rcf
 - Wash cell pellet once with 1x PBS
 - Resuspend pellet in 2-3 volumes of cell lysis buffer (at least 100 μ l)
 - Incubate on ice for 10 min
 - Centrifuge for 5 min at 5,000 rpm at 4°C, discard supernatant (cytoplasmic extract)
 - Resuspend nuclear pellet in 3-4 volumes of nuclei lysis buffer (at least 200 μ l, better more). If the pellet is difficult to resuspend use a syringe or sonify samples in a bioruptor (setting: 15 cycles, high, 30s on/off)
 - Incubate lysate for 30 min on a turning wheel at 4°C
 - Centrifuge for 10 min at 13,200 rpm and 4°C
 - Transfer supernatant to a fresh tube
 - Adjust lysate to 200 mM NaCl
 - Measure protein concentration (Bradford assay), do not use more than 2-3 mg of protein/experiment. If the protein is overexpressed 1 mg of lysate is enough
- *Capture of proteins by bead-bound RNA*
 - Add equal amounts of lysate to each aliquot of beads
 - Incubate on a turning wheel at 4 °C for 3-4 h
 - Wash beads 3x with nuclei lysis buffer
 - After the final wash add 20 μ l of 1x Laemmli buffer and heat samples for 10 min at to 75 °C.
 - Load samples onto a SDS PAA gel and detect proteins by Western blotting